LBM™
TSB Salt

Product Insert & How to Use Guide

PI09B rev.02 ; 2012.09

INTENDED USE

Copan TSB salt broth is a Triptic Soy Broth (TSB) with 2,5 % sodium chloride (NaCl). It is an enrichment medium able to support the growth of *S. aureus* spp. Microorganisms other than *S. aureus* may be supported or inhibited or partially inhibited.

Copan TSB salt broth is used to inoculate clinical specimens and, after appropriate incubation, is subcultured on selective agar plates.

SUMMARY

Copan TSB salt broth consists of Triptic Soy Broth (TSB) with 2,5% NaCl. Copan TSB salt is an enrichment medium used to support the growth of *S. aureus* spp., **Copan TSB salt broth** is used to inoculate clinical specimens and, after appropriate incubation, is subcultured on selective agar plates. Specific codes of Copan TSB salt broth can be used with Copan Walk-Away Specimen Processor (WASP), an automated system for planting and streaking of specimens onto culture media plates.

Copan TSB salt broth is available in different configurations: in bulk packs comprising multiple screw cap tubes filled with 2 ml of solution or as kits comprising of one screw cap tube containing 2 ml of solution plus three sterile regular size flocked swabs in peel pouch. TSB salt broth kit is specifically indicated for multi-site-sampling from the same patient.

REAGENTS**TSB salt broth solution components (per l):**

TSB powder : 30 g/l

Sodium Chloride powder : 25 g/l

PRECAUTIONS

1. For *in vitro* diagnostic use.
2.  This product is for single use only; reuse may cause a risk of infection and/or inaccurate results.
3. Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified personnel.
4. Work under a biological safety cabinet and wear gloves
5. All specimens and materials used to process them should be considered potentially infectious and handled in a manner which prevents infection of laboratory personnel. Sterilize all biohazard waste including specimen containers and media after their use.
6. Directions should be read and followed carefully.

STORAGE

This product is ready to use and no further preparation is necessary. The unopened bulk pack or individual kit pack can be stored at 5 - 25°C until used or until the expiration date. Do not overheat. Do not incubate or freeze prior to use. Improper storage will result in loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box.

PRODUCT DETERIORATION

Do not use **Copan TSB salt broth** if: (1) the product shows visible marks of damage or contamination; (2) there is evidence of leakage; (4) the expiration date has passed; (5) there are other signs of deterioration. (ie medium is turbid).

MATERIALS SUPPLIED

Catalog No.	Product Descriptions	Pack Size	Suitable for Automation
477CE.A	2 ml TSB with 2,5% NaCl medium in 12x80 mm polypropylene screw-cap tube with internal conical shape	50 units per shelf pack 6X50 units per box	YES
477CE02.A	2 ml TSB with 2,5% NaCl medium in 12x80 mm polypropylene screw-cap tube with internal conical shape plus 1 regular flocked swab and 1 regular flocked swabs pink shaft in peel pouch	50 units per shelf pack 6X50 units per box	YES
477CE02	2 ml TSB with 2,5% NaCl medium in 12x80 mm polypropylene screw-cap tube with internal conical shape plus 1 regular flocked swab white shaft and 1 regular flocked swabs pink shaft in peel pouch	50 units per shelf pack 6X50 units per box	NO
477CE03.A	2 ml TSB with 2,5% NaCl medium in 12x80 mm polypropylene screw-cap tube with internal conical shape plus 1 regular flocked swab white shaft and 2 regular flocked swabs pink shaft in peel pouch	50 units per shelf pack 6X50 units per box	YES
477CE03	2 ml TSB with 2,5% NaCl medium in 12x80 mm polypropylene screw-cap tube with internal conical shape plus 1 regular flocked swab with shaft and 2 regular flocked swabs pink shaft in peel pouch	50 units per shelf pack 6X50 units per box	NO

MATERIALS REQUIRED BUT NOT SUPPLIED

Appropriate materials for the cultivation and isolation of bacteria. Refer to laboratory reference manuals for recommended protocols for culture and identification techniques.

SPECIMEN COLLECTION**INSTRUCTIONS FOR USE****TSB SALT BROTH IN BULK****Direct specimen collection**

1. Take the tubes of TSB salt broth and unscrew the cap.
2. Transfer sample into the opened tube. Swab specimen should be broken or cut into the medium tube. For liquid specimens use a sterile loop to transfer the specimen into the medium tube.
3. Re-cap the tube of TSB salt broth and Vortex the tube for 5-10 seconds at 2000/2500 rpm in order to mix tube content.

The Copan ESwab system can be used for the specimen collection, refer to the PI of the ESwab product for instruction to use.

TSB salt broth kit :

1. Open the paper film pouch containing the tube of TSB salt broth and the paper film pouch containing the regular flocked swabs.
 2. Open the paper film pouch containing the regular flocked swabs.
 3. Take pink swab to collect first specimen (i.e. nose, throat or perineum sample) if you need to take two or more than two samples from the same patient.
- WARNING:** Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft. Swab shape may have breaking point and diameter change points that could result less resistant to pressure or bending. Do not use excessive force, pressure or bending when collecting swab sample in order to avoid accidental breakage. Unscrew the TSB salt broth tube.
4. Insert the swab into the TSB salt broth tube all the way to reach the bottom of the tube. Dip and gently stir the swab for 5 seconds.
 5. Lift up the swab from the liquid medium and swirl the swab against the tube walls 5 times to allow release of the sample from the flocked fibre. Remove the swab from the tube and recap.
 6. Discard the pink swab into Biohazard Waste.
 7. Repeat all previous step (3 to 6) and use a second pink flocked swab (if provided with the kit and) if you need to collect the second specimen from same patient (ie nose, throat or perineum sample). If not proceed to step 9.
- WARNING:** Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft. Swab shape may have breaking point and diameter change points that could result less resistant to pressure or bending. Do not use excessive force, pressure or bending when collecting swab sample in order to avoid accidental breakage.
8. Take white swab to collect the last specimen (ie nose, throat or perineum sample) and then break the swab at the moulded breaking point. Deposit the swab into the tube.
- WARNING:** Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft. Swab shape may have breaking point and diameter change points that could result less resistant to pressure or bending. Do not use excessive force, pressure or bending when collecting swab sample in order to avoid accidental breakage. Recap the tube with white swab inside, write patient name and send tube to the laboratory.
9. Vortex the tube for 5-10 seconds at 2000/2500 rpm in order to mix tube content.

SPECIMEN PROCESSING IN LABORATORY**PROCEDURE MANUAL OPERATION****TSB salt broth collected specimen**

1. Vortex the TSB salt inoculated broth for 10sec or according to the standard laboratory procedure

ESwab collected specimen

1. Take the tubes of TSB salt broth and unscrew the cap
2. Vortex the ESwab specimen tube for 10sec
3. Unscrew the cap and transfer the swab from the ESwab tube to the TSB salt tube by using sterile forceps; alternatively inoculate the TSB salt broth with 30µL of Amies inoculated media by a micropipette.
4. Re-cap the ESwab tube
5. Re-cap the tube of TSB salt broth
6. Vortex the inoculate TSB salt tube for 10sec

1. Incubate inoculated TSB salt broth tubes at 35 ± 2 °C .
2. Examine the TSB salt broth tube for growth after 18-24 hours.
3. Aseptically remove aliquots of the TSB salt broth and inoculate onto appropriate bacteriology selective medium.

NOTE: all the swabs could be broken into the tube. This procedure cannot be applied with 477CE03.A to be used in WASP. Moreover attention must be paid successive during tube decapping in order to avoid that the three swabs are contemporary extracted.

PROCEDURE AUTOMATIC OPERATION

The TSB Salt broth is manufactured in a container suitable to be processed by automatic system for planting and streaking (WASP). By using the WASP is possible to inoculate specimens collected by the ESwab system into the broth and to subculture the enriched broth sample.

BROTH INOCULATION BY THE WASP

1. Process the specimen collected in ESwab by a specific broth inoculation protocol by using the 30µL loop. Refer to the WASP User Manual for further information.
2. Incubate the broth at 35°C ±2°C for 18-24 hours.

SUBCULTURE OF THE ENRICHED BROTH SAMPLE BY THE WASP

1. Process the enriched TSB salt sample by a specific protocol for the streaking of an appropriate selective medium plate. The use of a 1 or 10µL loop and a multiple quadrant streaking pattern (3-4 Quadrants type) is suggested in order to improve the colonies isolation. It is recommended the laboratory to validate the streaking protocol. Refer to the WASP User Manual for further information.
2. Incubate the plates according to the standard laboratory procedure.

LIMITATIONS

1. In the laboratory, wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens.
2. Condition, timing, and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection.
3. Performance testing with Copan TSB salt broth was conducted using laboratory strains spiked into the TSB salt broth tube and not using human specimens.
4. Proper specimen collection from the patient is extremely critical for successful isolation and identification of infectious organisms. For specific guidance regarding specimen collection procedures, consult published reference manuals. Specimens should be collected as soon as possible after the clinical onset of disease. Highest bacterial titres are present during the acute illness.
5. Microorganisms other than S aureus may be supported or inhibited or partially inhibited.

WARNINGS

1. Do not re-pack.
2. Not suitable for any other application than intended use.
3. The use of this product in association with any diagnostic assay or with any diagnostic instrumentation should be validated by the user before using.
4. Do not use if the product is visibly damaged
5. Do not ingest the medium.
6. Directions for use must be followed carefully. The manufacturer cannot be held responsible for any unauthorized or unqualified use of the product.
7. Product to be handled by trained personnel only.
8. It must be assumed that all specimens contain infectious micro-organisms; therefore all specimens must be handled with appropriate precautions. After use, tubes must be disposed of according to laboratory regulations for infectious waste.
9.  Copan TSB salt broth is for single use only; reuse may cause a risk of infection and/or inaccurate results.
10. Copan TSB salt broth is for in-vitro diagnostics use only and is in no way intended for a curative or prophylactic purposes.

WASTE DISPOSAL

Unused reagents may be considered as non hazardous waste and disposed of accordingly.
Dispose of used reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products. It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

RESULTS:

Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the laboratory.

PERFORMANCE CHARACTERISTICS AND QUALITY CONTROL PROCEDURE:**Procedure for testing:**

1. Starting from a fresh culture plate of tester organism prepare 0.5 Mc Farland inoculum of the tester organism in PBS
2. prepare appropriate dilutions in PBS of original 0.5 McF inoculum in order to obtain a diluted inoculum containing from 30 to 300 CFU/100 µl of TSB
3. with a micropipette spike 100 µl of chosen dilution into TSB salt broth tube
4. vortex inoculated tube for 10 seconds at 2500-3000 RPM
5. spread 100 µl of inoculated TSB salt broth tube onto MSA for the zero time count
6. incubate the inoculated plate at 35+/-2°C for 18-24 hours
7. incubate the inoculated TSB salt broth tube at 35+/-2°C for 18-24 hours
8. vortex inoculated tube for 10 seconds at 2500-3000 RPM
9. spread 100 µl of inoculated TSB salt broth tube onto MSA (Mannitol Salt Agar)
10. incubate the inoculated plate at 35+/-2°C for 18-24 hours

ACCEPTABILITY LIMITS : GROWTH.**PERFORMANCE TEST RESULTS:**

STRAIN [‡]	ZERO TIME COUNT; CFU/PLATE	CFU COUNT AFTER 24 HOURS at 35°C±2°C; CFU/PLATE
MRSA ATCC 43300	126	CONFLUENT GRWOTH
S.aureus ATCC 6538	112	CONFLUENT GRWOTH

[‡] Performance testing with Copan TSB salt broth was conducted using laboratory strains. Performance testing was not conducted using human specimens.

BIBLIOGRAPHY

- Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger and W.C. Winn, Jr. 1992. Color Atlas and Textbook of Diagnostic Microbiology. 4th ed. J.B. Lippincott Co. Philadelphia, PA.
- Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology. 7th ed. ASM, Washington, D.C..
- National Committee for Clinical Laboratory Standards (NCCLS). 2003. Quality Control of Microbiological Transport Systems. Approved Standard M40-A.
- 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
- Miller, J. M. 1999. A Guide to Specimen Management in Clinical Microbiology, 2nd ed. ASM, Washington, DC.
- Isenberg, H. D., 2004. Clinical Microbiology Procedures Handbook, 2nd ed. ASM, Washington, DC.
- Isenberg, H.D., 1998. Essential Procedures for Clinical Microbiology. Chapter 14.12, Page 787. Packaging and Shipping Infectious Substances.
- National Committee for Clinical Laboratory Standards (NCCLS). 1994. Procedures for Handling and Transport of Diagnostic Specimens and Etiologic Agents. Approved Standard H5-A3.

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BIBLIOGRAFIA

- Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger and W.C. Winn, Jr. 1992. Color Atlas and Textbook of Diagnostic Microbiology. 4th ed. J.B. Lippincott Co. Philadelphia, PA.
- Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology. 7th ed. ASM, Washington, D.C..
- National Committee for Clinical Laboratory Standards (NCCLS). 2003. Quality Control of Microbiological Transport Systems. Approved Standard M40-A.
- 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
- Miller, J. M. 1999. A Guide to Specimen Management in Clinical Microbiology, 2nd ed. ASM, Washington, DC.
- Isenberg, H. D., 2004. Clinical Microbiology Procedures Handbook, 2nd ed. ASM, Washington, DC.
- Isenberg, H.D., 1998. Essential Procedures for Clinical Microbiology, Chapter 14.12, Page 787. Packaging and Shipping Infectious Substances.
- National Committee for Clinical Laboratory Standards (NCCLS). 1994. Procedures for Handling and Transport of Diagnostic Specimens and Etiologic Agents. Approved Standard H5-A3.

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BIBLIOGRAPHIE

- Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger and W.C. Winn, Jr. 1992. Color Atlas and Textbook of Diagnostic Microbiology. 4th ed. J.B. Lippincott Co. Philadelphia, PA.
- Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology. 7th ed. ASM, Washington, D.C.
- National Committee for Clinical Laboratory Standards (NCCLS). 2003. Quality Control of Microbiological Transport Systems. Approved Standard M40-A.
- 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
- Miller, J. M. 1999. A Guide to Specimen Management in Clinical Microbiology, 2nd ed. ASM, Washington, DC.
- Isenberg, H. D., 2004. Clinical Microbiology Procedures Handbook, 2nd ed. ASM, Washington, DC.
- Isenberg, H.D., 1998. Essential Procedures for Clinical Microbiology. Chapter 14.12, Page 787. Packaging and Shipping Infectious Substances.
- National Committee for Clinical Laboratory Standards (NCCLS). 1994. Procedures for Handling and Transport of Diagnostic Specimens and Etiologic Agents. Approved Standard H5-A3.

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
Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger and W.C. Winn, Jr. 1992. Color Atlas and Textbook of Diagnostic Microbiology. 4th ed. J.B. Lippincott Co. Philadelphia, PA.
 Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology. 7th ed. ASM, Washington, D.C..
 National Committee for Clinical Laboratory Standards (NCCLS). 2003. Quality Control of Microbiological Transport Systems. Approved Standard M40-A.
 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
 Miller, J. M. 1999. A Guide to Specimen Management in Clinical Microbiology, 2nd ed. ASM, Washington, DC.
 Isenberg, H. D., 2004. Clinical Microbiology Procedures Handbook, 2nd ed. ASM, Washington, DC.
 Isenberg, H.D., 1998. Essential Procedures for Clinical Microbiology. Capítulo 14.12, Página 787. Packaging and Shipping Infectious Substances.
 National Committee for Clinical Laboratory Standards (NCCLS). 1994. Procedures for Handling and Transport of Diagnostic Specimens and Etiologic Agents. Approved Standard H5-A3.

Index of Symbols/ Table des Symboles/Symbole/Tabla de Símbolos/Tabella dei Simboli/Quadro de Símbolos

Symbol/Symbole/Símboło	Meaning/Signification/Bedeutung/Significado/Significato/Significado
	Manufacturer/Fabricant/Hersteller/Fabricante/Fabbricante/Fabricante
	In Vitro Diagnostic
	Do no reuse/Ne pas réutiliser/Nicht zur Wiederverwendung/No reutilizar/Non riutilizzare/Não voltar a usar
	Catalogue number/Référence du catalogue/Bestellnummer/Número de catálogo/Numero di catalogo/Referência do catálogo
	Temperature limitation/Limites de temperature/Temperaturbegrenzung/Límites de temperatura/Limiti di temperatura/Límites de temperatura
	Use by/Utiliser jusque/Verwendbar bis/Fecha de caducidad/Utilizzare entro/Prazo de validade
	Consult Instructions for Use/Consulter les instructions d'utilisation/Gebrauchsanweisung beachten/Consultar las instrucciones de uso/Consultare le istruzioni per l'uso/Consultar as instruções de utilização
	Peel/Décoller/Abziehen/Desprender/Strappare per aprire/Destacável
	Batch code (Lot)/Code de lot (Lot)/Chargencode (Chagenbezeichnung)/Código de lote (Lote)/Codice del lotto (partita)/ Código do lote (Lote)
	Contains sufficient for <n> tests/Contenu suffisant pour <n> tests/Ausreichend für <n> Tests/Contenido suficiente para <n> pruebas/Contenuto sufficiente per <n> test/Contém suficiente para <n> testes
	<p>Copan Italia Spa Via Perotti 10 25125 Brescia Italia</p> <p>Teléfono: +39 030 2687211 Fax: +39 030 2687250 E-mail: info@copanitalia.com Página web: www.copanitalia.com</p> 